



Science.  
Applied to Life.™



## 3M™ Imprint™ 4 Bite VPS Bite Registration Material

# Fast and accurate bite registrations

With its great flow properties and no-slump consistency, Imprint 4 Bite material does not lock into undercuts or interdental areas. Material hardness and low flexibility after setting allows for easy cutting and trimming as well as accurate positioning of stone models.

### **Saves time with a fast intra-oral setting.**

Imprint 4 Bite material offers a convenient working time and a short setting time of only 60 seconds, keeping procedures moving quickly while enhancing patient comfort.



# Indications

3M™ Imprint™ 4 Bite Registration material is a universal VPS bite registration material. The material can also be scanned for occlusal recording and registration with CAD/CAM systems.



Scan of bite registration made with Imprint 4 Bite material.

	Mixing and application	Setting time from start of mixing*
<b>3M™ Imprint™ 4 Bite VPS Bite Registration Material</b>	0:20	1:00


\*The stated times are valid at 23°C/74°F. Times depend on ambient temperature.

## Get more confidence in every step.

### 1 Retract



### 2 Syringe and Mix



### 3 Impress



### Preliminary Impression



### Bite Registration



## Ordering information

Item No.	Product Information
71529	<b>3M™ Imprint™ 4 Bite Refill</b> 2 cartridges – 50ml each, 5 Garant™ Mixing Tips green, 5 Garant™ Intra-oral Tips Bite
71460	<b>3M™ Garant™ Intra-oral Tips Bite</b> Includes 50 Intra-oral Tips Bite



3M™ Imprint™ 4 Bite Refill (71529)



3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249  
Web 3M.com/dental

3M Canada  
Post Office Box 5757  
London, Ontario N6A 4T1  
Canada

Phone 1-888-363-3685

Customer Care Center: 1-800-634-2249  
3M.com/greatimpressions

3M, Garant and Imprint are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Please recycle. Printed in U.S.A.  
© 3M 2016. All rights reserved. 70-2013-0627-4 (Rev. 8/2016)