

## AH Temp™ Root Canal Dressing



## AH Plus Jet® Root Canal Sealer



# The perfect match!

AH Temp™ combines effective disinfection with easy handling

## Goal of Endodontic Treatment: Effective Disinfection of Root Canal


AH Temp™ is a highly alkaline paste with a pH level of 12 based on calcium hydroxide which is known for its antibacterial effect.

### Reduced Desiccation

- Less blocking of cannulas during application
- Less blocking of syringes over time
- Less waste of desiccated syringes


#### Stable Consistency during Chairside Use<sup>1</sup>

[hrs.]




Once the syringe cap is removed, AH Temp™ keeps its pasty consistency for 6 hours – no time pressure during the treatment.

In the daily routine, AH Temp™ is usable for 3 months



### Excellent Radiopacity enables reliable X-ray Control<sup>1</sup>

[mm/mm Al]



<sup>1</sup> Data on file.

<sup>2</sup> Not a registered trademark of Dentsply Sirona, Inc.

# AH Temp™ – Direct and Precise Application into the Root Canal



“The application of AH Temp™ is brilliant – it is the best root canal dressing available at the moment. In addition, it is supplied with an excellent cannula!”

Dr. Matthias J. Roggendorf,  
University of Marburg, Germany



AH Temp™ Cannula




Calxyl² Cannula

## Proven Compatibility with AH Plus Jet®

### Consistent High Bond Strength of AH Plus Jet® on Dentin<sup>1</sup>

[MPa]



#### The perfect match!

AH Temp™ is compatible with AH Plus Jet® – the ideal combination of root canal dressing and sealer.

Proven consistent performance of AH Plus Jet® after the removal of AH Temp™ from the root canal.



## Handling Benefits of AH Temp™ and AH Plus Jet®

- AH Temp™ – Calcium Hydroxide-Based Root Canal Dressing with ultra-thin, long and pre-angulated AH Temp™ Cannulas
- Automatic mixing for convenient, easy use.
- It's the safe and convenient way to achieve precise placement, even in difficult anatomic conditions.



<sup>1</sup> M. Roggendorf, 2012, proven for application of AH Temp™ and AH Plus Jet® Sealer, data on file.

<sup>2</sup> Not a registered trademark of Dentsply Sirona, Inc.


<sup>3</sup> Data on file.

# Ordering Information

## AH Temp™ Kit

4 x AH Temp™ Syringes, 0.75 ml each,  
20 x AH Temp™ Cannulas,  
1x Finger Grip,  
1x Illustrated Technique Guide,  
1x Instructions For Use

606.20.123



## AH Temp™ Cannula Refill

20 x AH Temp™ Cannulas

606.20.124



## AH Plus Jet® Starter Kit

1 x AH Plus Jet® Mixing Syringe, 15 g,  
20 x Mixing Tips,  
20 x Intra Oral Tips  
1 x Mixing Pad,  
1 x Illustrated Technique Guide,  
1 x Organizer Tray

606.20.115



## AH Plus Jet® Refill

2 x AH Plus Jet® Mixing Syringe, 15 g

606.20.118



## MixingTips/Intra Oral Tips for AH Plus Jet®

40 x Mixing Tips with Intra Oral Tip

606.20.116



For more information visit [www.dentsplysirona.com](http://www.dentsplysirona.com).

Find us on  [www.facebook.com/dentsplysirona.endodontics](https://www.facebook.com/dentsplysirona.endodontics)

K79200094-05  
© Dentsply Sirona 07/2019